

FREE PROSTATE SPECIFIC ANTIGEN (f-PSA) ENZYME IMMUNOASSAY TEST KIT

IVD For in vitro diagnostic and professional use only

2°C  8°C
Store at 2 to 8°C.

 96 Tests

INTENDED USE

For the quantitative determination of f-PSA concentration in human serum.

INTRODUCTION

Prostate Specific Antigen PSA, is a 33KD serine protease, a single chain glycoprotein is an enzyme found in the blood produced exclusively by prostate cells. Normal levels of PSA in the blood are small amounts between 0-2.5 ng/ml. Higher than normal levels, greater than 2.5 ng/ml, can be caused by cancer or benign, non-cancerous conditions such as enlarged prostate, prostate inflammation, infection, or trauma. All elevated readings of PSA should be checked.

Current methods of screening men for prostate cancer utilize the detection of the major PSA-ACT form. Levels of 4.0 ng/ml or higher are strong indicators of the possibility of prostatic cancer. However, elevated serum PSA levels have also been attributed to benign prostatic hyperplasia and prostatitis, leading to a large percentage of false positive screening results. A potential solution to this problem involves the determination of free PSA levels.

PRINCIPLE OF THE TEST

The f-PSA ELISA test is based on the principle of a solid phase enzyme-linked immunosorbent assay. Anti-f-PSA monoclonal is coated on the microtiter wells and another anti-PSA monoclonal antibody labeled with horseradish peroxidase is used as the tracer. The f-PSA molecules present in the standard solution or serum are "sandwiched" between the two antibodies. Following the formation of the coated antibody-antigen-antibody-enzyme complex, the unbound antibody-enzyme tracers are removed by washing. The horseradish peroxidase activity bound in the wells is then assayed by a colorimetric reaction. The intensity of the color formed is proportional to the concentration of f-PSA present in the sample.

MATERIALS PROVIDED

- Anti-f-PSA coated microtiter plate with **96 wells**.
- Sample diluent (zero buffer), **12 ml transparent cap bottle with yellow**

color reagent.

- Reference standard set containing :
0 (silver cap vial) , 0.1 (golden cap vial), 0.5 (green cap vial) , 2.0 (blue cap vial) ,5.0 (purple cap vial) , and 10.0 (red cap vial) ng/ml ,(liquid ready to use) or lyophilized form, **1ml** each.
- Enzyme Conjugate Reagent, **22 ml red cap bottle.**
- TMB Reagent solution, **12 ml brown cap bottle.**
- Stop Solution, **12 ml white cap bottle.**
- Wash buffer concentrate (50X), **15ml green cap bottle.**
- Control set (**optional**).

Materials required but not provided:

- Vortex mixer .
- Precision pipettes:0.10,0.20 and 10.0ml.
- Disposable pipette tips.
- Graph paper.
- Distilled water.
- Absorbent paper or paper towel.
- A microtiter plate reader with a bandwidth of 10 nm or less and an optical density range of 0-2 OD or greater at 450 nm.

SPECIMEN COLLECTION AND PREPARATION

1. Blood should be drawn using standard venipuncture Techniques and the serum should be separated from the red blood cells as soon as possible. Avoid grossly hemolytic, Lipemic or turbid samples.
2. Plasma samples collected in tubes containing EDTA, Heparin, or oxalate may interfere with the test procedures and should be avoided.
3. Specimens should be capped and may be stored up to 48 hours at 2-8°C, prior to assaying. Specimens held for a longer time can be frozen at -20°C. Thawed samples must be mixed prior to testing.

STORAGE OF TEST KIT

1. Unopened test kits should be stored at 2-8°C upon receipt. The microtiter plate should be kept in a sealed bag with desiccants, to minimize exposure to damp air.
2. Opened test kits will remain stable until the expiration date, provided it is stored as described above.
3. A microtiter plate reader with a band width of 10 nm or less and an optical density (O.D) range of 0-2.5 or greater at 450 nm Wavelength is acceptable for use in absorbance measurement.

REAGENT PREPARATION

1. All reagents should be brought to room temperature (18-22°C) and mixed by gently inverting or swirling prior to use. Do not induce Foaming.
2. Dilute 1 volume of Wash Buffer (50x) with 49 volumes of Distilled water. For example, Dilute 15 mL of Wash Buffer (50x) into

distilled water to prepare 750 mL of washing buffer (1x). Mix well before use.

ASSAY PROCEDURE

1. Secure the desired number of coated wells in the holder.
2. Dispense 100 µl of standards, specimens, and controls into appropriate wells.
3. Dispense 100 µl of sample diluent into each well.
4. Thoroughly mix for 10 seconds. It is very important to have a complete mixing in this step.
5. Incubate at room 37°C for 60 minutes.
6. Remove the incubation mixture by emptying plate contents into a waste container.
7. Rinse and empty the microtiter wells 5 times withwashing buffer (1X).
8. Strike the wells sharply onto absorbent paper or paper towels to remove all residual water droplets.
9. Dispense 200 µl of Enzyme Conjugate Reagent into each well. Gently mix for 5 seconds.
10. Incubate at 37°C 60 minutes.
11. Remove the incubation mixture by emptying plate contents into a waste container.
12. Rinse and empty the microtiter wells 5 times with wash buffer (1X).
13. Strike the wells sharply onto absorbent paper to remove residual water droplets.
14. Dispense 100 µl of TMB Reagent into each well. Gently mix for 5 seconds.
15. Incubate at room temperature for 20 minutes.
16. Stop the reaction by adding 100 µl of Stop Solution to each well.
17. Gently mix for 30 seconds. It is important to make sure that all the blue color changes to yellow color completely.
18. Using a microliter plate reader, read the optical density at 450nm within 15 minutes.

WASH PROCEDURE

1. The wash procedure is critical. Insufficient washing will Result in a poor precision and falsely elevated absorbance readings.
2. It is recommended that no more than 32 wells be used for each assay run, if manual pipetting is used, since pipetting of all standards, specimens and controls should be Completed within 3 Minutes. A full plate of 96 wells may be used if automated pipetting is available.
3. Duplication of all standards and specimens, although not required, is recommended.

CALCULATION OF RESULTS

1. Calculate the average absorbance values (A_{450}) for each set of

reference standards, control, and samples.


- Construct a standard curve by plotting the mean absorbance obtained for each reference standard against its concentration in ng/ml on linear graph paper, with absorbance on the vertical (y) axis and concentration on the horizontal (x) axis.
- Use the mean absorbance value for each sample, determine the corresponding concentration of PSA in ng/ml from the standard curve.

EXAMPLE OF STANDARD CURVE

Results of a typical standard run with optical density readings at 450nm shown in the Y axis against f-PSA concentrations shown in the X axis. This standard curve is for the purpose of illustration only, and should not be used to calculate unknowns. Each user should obtain his or her own data and standard curve.

F-PSA (ng/ml)	Absorbance (450 nm)
0.0	0.006
0.1	0.032
0.5	0.155
2.0	0.597
5.0	1.302
10.0	2.361

O.D



EXPECTED VALUES AND SENSITIVITY

Clinical studies have shown that higher total PSA levels and lower percentages of free PSA are associated with higher risks of prostate cancer. of free PSA are associated with higher risks of prostate cancer. In one study, 52 individuals with benign prostate hyperplasia (BPH) and 77 individuals with prostate cancer were analyzed. The f-PSA/t-PSA ratio were summarized in the table below.

	f-PSA/t-PSA Ratio			
	Median	Min	Max	Mean (95% confidence interval)
BPH	0.18	0.04	0.42	0.19 (0.17-0.21)
Prostate Cancer	0.09	0.02	0.53	0.12 (0.10-0.14)

The relationship between f-PSA/t-PSA ratio and risk of prostate cancer is also age related. When total PSA is in the range of 4.0-10.0 ng/mL, a f-PSA/t-PSA ratio ≤ 0.10 indicates 49% to 65% risk of prostate cancer depending on age ; a f-PSA/t-PSA ratio > 0.25 indicates a 9% to 16% risk of prostate cancer, depending on age.

f-PSA/t-PSA ratio	Percent probability of prostate cancer		
	50-59 years	60-69 years	=70 years
≤ 0.10	49.2	57.5	64.5
0.11-0.18	26.9	33.9	40.8
0.19-0.25	18.3	23.9	29.7
> 0.25	9.1	12.2	15.8

Multiple factors such as population, age, specificity of test method may affect interpretation of f-PSA and t-PSA values. The seranges should be used as guidelines only. Each laboratory should establish its own reference values.

Sensitivity

The minimum detectable concentration of f-PSA in this assay is estimated to be 0.05 ng/ml.

Cross Reactivity

Antigens	Concentration	%Cross-react
PSA-ACT	500 ng/mL	0.2
AFP	10,000 ng/ mL	0
CEA	5,000 ng/ mL	0
CA 125	1,000 U/ mL	0
CA 15-3	1,000 U/ mL	0
CA 19-9	1,000 U/ mL	0
-HCG	1,000 ng/ mL	0
β -HCG	1,000 ng/ mL	0
HCG	50,000 mIU/ mL	0

Precision

Intra- Assay			
	Replicates	S.D.	%CV
level	20	0.005	13.1
level	20	0.011	4.4
level	20	0.128	3.2

Linearity

Two patient sera were serially diluted with 0 ng/ml standard. The average recovery was 108.7%.

Sample A

Dilution	Expected	Observed	% Recov
undiluted	8.691	8.691	100
2X	4.346	4.455	102.5
4X	2.173	2.290	105.4
8X	1.086	1.258	115.8
16X	0.543	0.617	113.6
32X	0.272	0.294	108.1
64X	0.136	0.147	108.1

Average Recovery: 107.6%

Sample B

Dilution	Expected	Observed	% Recovery
undiluted	7.015	7.015	100
2X	3.508	3.516	100.2
4X	1.754	1.834	104.6
8X	0.877	0.970	110.6
16X	0.438	0.509	116.2
32X	0.219	0.249	113.7
64X	0.110	0.136	123.6

Average Recovery: 109.8%

Recovery

Equal parts of diluted patient sera were mixed to test for interference by unknown materials, such as drugs or hormones, in the assay. Concentrations of Free PSA were determined before (original and added) and after (observed). The average recovery was 99.2%.

Sample 1

Sample	Orig.Conc	Added	Expected	Obsearve d	Recovery
A	9.802	0.141	4.972	4.540	91.3
B	9.802	0.281	5.042	4.638	92.0
C	4.699	2.168	3.434	3.342	97.3
D	2.168	1.170	1.669	1.661	99.5
E	0.563	0.281	0.422	0.423	100.2
F	4.699	1.125	2.912	2.727	93.6
G	0.563	0.141	0.352	0.399	113.4

Average Recovery: 98.2%

Sample 2

Sample	Orig.Conc	Added	Expected	Obsearve d	Recovery
A	6.825	0.098	3.462	3.301	95.3
B	6.825	0.195	3.510	3.156	89.9
C	3.342	1.563	2.453	2.503	102.0
D	1.706	0.781	1.244	1.205	96.9
E	0.391	0.195	0.293	0.327	111.6
F	3.342	0.781	2.062	1.964	95.2

G	0.391	0.098	0.245	0.271	110.6
---	-------	-------	-------	-------	-------

Average Recovery: 100.2%

Limitations of the Procedure

There are some limitations of the assay:

1. As with all diagnostic tests, a definite clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
2. Studies have implicated possible interference in immunoassay results in some patients with known rheumatoid factor and antinuclear antibodies. Serum samples from patients who have received infusions containing mouse monoclonal antibodies for diagnostic or therapeutic purposes, may contain antibody to mouse protein (HAMA). Although we have added some agents to avoid the interferences, we cannot guarantee it will eliminate all the effects of that.

REFERENCES

1. Christensson A, Bjork T, Nilsson O, et al. Serum Prostate Specific Antigen Complexed to α 1-Antichymotrypsin As An Indicator of Prostate Cancer. J. of Urol. 150:100-105; 1993.
2. Lilja H, Christensson A, Dahlen U, et al. Prostate-specific antigen in serum occurs predominantly in complex with α 1-antichymotrypsin. Clin Chem. 1991;37:1618-1625.
3. Stenman U-H, Leinonen J, Alfthan H, Rannikko S, Tuhkanen K, Alfthan O.A complex between prostate-specific antigen and α 1-antichymotrypsin is the major form of prostate-specific antigen in serum of patients with prostate cancer: assay of the complex improves clinical sensitivity for cancer. Cancer Res. 1991;51:222-226.
4. Catalona WJ, Smith DS, Ratliff TL, Basler JW. Detection of organ-confined prostate cancer is increased through prostate-specific antigen-based screening. JAMA. 1993;270:948-954.
5. Stamey TA, Yang N Hay AR, McNeal JE, Freiha, FS, Redwine E. Prostate-specific antigen as a serum marker for adenocarcinoma of the prostate. N Engl J Med. 1987;317:909-916.
6. Catalona et al. Percentage of Free Serum PSA and Prostate Cancer Detection. JAMA. 1995;274, No. 15: 1214-1220.
7. Catalona WJ, Smith DS, Ratliff TL, et al. Measurement of prostate-specific antigen in serum as a screening test for prostate cancer. N Engl J Med. 1991;324:1156-1161. Erratum: N Engl J Med. 1991;325:1324.
8. Carter HB, Pearson JD, Metter J, et al. Longitudinal evaluation of prostate-specific antigen levels in men with and without prostate disease. JAMA. 1992;267:2215-2220.
9. Smith DS, Catalona WJ. Rate of change in serum prostate-specific antigen levels as a method for prostate cancer detection. J Urol. 1994;152:1163-1167.
10. Benson MC, Whang IS, Pantuck A, et al. Prostate specific antigen density: a means of distinguishing benign prostatic hypertrophy and prostate cancer. J Urol. 1992; 147:815-816.

11. Catalona WJ, Hudson MA, Scardino PT, et al. Selection of optimal prostate specific antigen cut-offs for early detection of prostate cancer: receiver operating characteristic curves. J Urol. 1994;152:2037-2042.
12. Smith DS, Catalona WJ. The nature of prostate cancer detected through prostate specific antigen based screening. J Urol. 1994;152:1732-1736.
13. Oesterling JE, Jacobsen SJ, Chute CG, et al: Serum prostate-specific antigen in a community-based population of healthy men. JAMA 1993 Aug 18;270:860-864.
14. Catalona WJ, Smith DS, Wolfert RL, et al: Evaluation of percentage of free serum prostate-specific antigen to improve specificity of prostate cancer screening. JAMA 1995;274(15);214-1220.
15. Jacobsen SJ, Bergstralh EJ, Guess HA, et al: Predictive properties of serum prostate-specific antigen testing in a community-based setting. Arch Intern Med 1996;156:2462-2468.
16. Oesterling JE, Jacobsen SJ, Klee GG, et al: Free, complexed and total serum prostate specific antigen: the establishment of appropriate reference ranges for their concentrations and ratios. J Urol 1995;154:1090-1095.
17. Dworschack RT, Thiel RP, Picolli SP: Clinical evaluation of the free/total PSA ratio generated with the Elecsys total and free PSA assays on the Elecsys 1010 and 2010 systems. Clin Chem 2001 Jun; 47 (6 Suppl):A149 (Abstract).



ATLAS Medical

Ludwig-Erhard Ring 3

15827 Blankenfelde-Mahlow

Germany

Tel: +49 - 33708 – 3550 30

Email: Info@atlas-medical.com

PPI1660A01

Rev A (02.09.2019)

	Catalogue Number		Temperature limit
	In Vitro diagnostic medical device		Caution
	Contains sufficient for <n> tests and Relative size		Consult instructions for use (IFU)
	Batch code		Manufacturer
	Fragile, handle with care		Use-by date
	Manufacturer fax number		Do not use if package is damaged
	Manufacturer telephone number		Date of Manufacture
	Keep away from sunlight		Keep dry